


# Operating instructions

## Density determination set for precision balance KERN EMB 200-3V

### KERN YDB-01

Version 1.0  
04/2012  
GB



YDB-01-BA-e-1210



# KERN YDB-01

Version 1.0 04/2012

Operating instructions

Density determination set for precision balance KERN EMB 200-3V


## 1 Scope of delivery

- ⇒ Check packaging and density determination set immediately when unpacking for possible visible damage.
- ⇒ Make sure that all parts are completely present.

- 1 Weighing plate "density set"



- 2 Platform



- 3 Beaker



- 4 Immersion basket for descending solid matter (density > 1 g/cm<sup>3</sup>)



- 5 Immersion basket for floating solid matter (density < 1 g/cm<sup>3</sup>)






- 6 Sinker  
20 g stainless steel weight



- 7 Thermometer





## 2 Dimension [mm]

| Immersion basket for descending solid matter (density $> 1 \text{ g/cm}^3$ )   | Immersion basket for floating solid matter (density $< 1 \text{ g/cm}^3$ )  |
|--|---|
|  <p>A technical drawing of an immersion basket for descending solid matter. It consists of a top flared rim, a central vertical rod with a small hole near the top, and a bottom flared rim. A vertical dimension line on the right indicates a total height of (115) mm.</p>    |  <p>A technical drawing of an immersion basket for floating solid matter. It features a top flared rim, a central vertical rod with a small hole near the top, and a wider, flat bottom flared rim. A vertical dimension line on the right indicates a total height of (101) mm.</p> |
| <p style="text-align: center;">Sample dish</p>  <p>A technical drawing of a sample dish. It is a shallow, wide, flared bowl with a central stem. The top diameter is labeled as <math>\phi 35</math> and the bottom diameter of the stem is labeled as <math>\phi 9</math>.</p> |   |

Beaker



Platform



### 3 Installation with KERN EMB 200-3V



Fig.1: KERN EMB 200-3V with installed density sets KERN YDB-01 (numbering, see chap. 1)

- ⇒ Disconnect scale from power supply.
- ⇒ Remove the standard weighing plate and replace it with the density set.
- ⇒ Place the platform for glass containers in a way that it does not touch the weighing plate.
- ⇒ Place beaker in the centre of the platform. Make sure that it has no contact with the frame.
- ⇒ Hang the immersion basket on the rack. Make sure that it is centred in the recess.
- ⇒ Pour the liquid into the glass beaker. Filling height should be approx.  $\frac{3}{4}$  of the capacity. Immerse thermometer
- ⇒ Heat the liquid, the instruments or the sinker until the temperature is constant. Observe the warm-up time of the balance.



For further information, as well as the procedure of density determination please see the operating instructions added to the density balance KERN EMB 200-3V or on the KERN homepage ([www.kern-sohn.com](http://www.kern-sohn.com)).